

DOSI-FUSER[®] *HIGH VOLUME*

400
500
600



Safety and Trust

The widest range of sizes and flow rates

The same Main Product Features than DOSI-FUSER® Range

- **Unique DOSI-FUSER®- DESIGNED CAPILLARY**
- **Specially designed guided balloon with Scale Indicator**
- **Product range with easy identification by coloured labels**
- **Easy-to-fill system**


Features, Advantages & Benefits

RESERVOIR

- **Rigid plastic housing and UV-light protection**
- **Polyisoprene balloon reservoir**
- **Independent entry and exit ports**

INFUSION LINE

- **Infusion line to patient:**
kink-resistant with removable clamp
- **Safety in-line filters:** 1.2 µm particle filter with ≤0.03 µm air-removal filter



DOSI-FUSER® HIGH VOLUME	400 ml	500 ml	600 ml
2 days	2days 400ml 8.3ml/h 25915-400D2	2days 500ml 10.4ml/h 25915-500D2	2days 600ml 12.5ml/h 25915-600D2
1 day			1day 600ml 25ml/h 25915-600D1
12 hours		12hours 500ml 41.6ml/h 25915-500H12	12hours 600ml 50ml/h 25915-600H12
4 hours	4hours 400ml 100ml/h 25915-400H4		
2 hours		2hours 500ml 250ml/h 25915-500H2	

This product is certified by the Notified Body, UL International (UK) Ltd (0843), in accordance with Council Directive 93/42/EEC of Medical Devices

