


VeinCLEAR™ System Catheter



Heating Element Length



VeinCLEAR™ System offers

- Outpatient procedure
- Minimum complications
- High occlusion rate
- Short recovery time

Features

- 5Fr and 7Fr catheters
- Heating element length from 1cm to 7cm
- Available in Alloy52, Stainless Steel, LED and Guidewire options

VeinCLEAR™ System

Endovenous Radiofrequency Ablation System



RF Medical Co., Ltd.

#503, 506, 254, Beotkkot-ro, Geumcheon-gu, Seoul, Korea

Tel +(82) 2-2108-4200 Fax +(82) 2-2113-7999 E-mail overseas@rfa.co.kr Homepage www.rfa.co.kr

RFCT-003(E) Rev. 7 (2019.09.30)



01 VeinCLEAR™ System

Dedicated for Endovenous Thermal Ablation

+ Generator V-700

- LCD touch screen
- Compatible with all STS and Alloy52 catheters
- Advanced automatic self-test system
- Safe temperature control algorithm
- PC compatibility for observation and data download
- Multi-language support with audio guidance
- Target temperature and ablation cycle time can be optionally adjusted
- Pre-set ablation algorithm available according to vein diameter



Specifications

Product Code	V-700
Protection Class	I
Protection Type	BF
Input Power Voltage	100 - 240V (AC)
Input Power Frequency	50/60 Hz
Maximum Input Power	150 VA
Measuring Temperature	0°C - 200°C
Alarm Sound	More than 65 dB

* M-3004 and V-1000 generators are compatible with Alloy52 catheters only.

02 Tumescentor™ System

Tumescent Infiltration for Endovenous Thermal Ablation

+ Features

- Operation status display monitor
- Injection speed control
- Dedicated IV pole included
- Tumescentor™ accessories are compatible with:
 - RF Medical RFP-300 pump (hand piece option only)
 - Nouvag DP-30 pump (simple tubing & hand piece option)




Specifications


Product Code	RFP-300
Protection Class	I
Protection Type	BF
Input Power Voltage	100 - 240V (AC)
Input Power Frequency	50/60 Hz
Maximum Input Power	100 VA

03 VeinCLEAR™ Catheters

Dedicated for Endovenous Thermal ablation



+ Options



Specifications

Product Code	Compatible Introducer	Heating Element	Total Length	Heating Element Length	Compatible Guidewire	
VVT 70FS	5Fr	Stainless Steel	70cm			
VVT 100FS			100cm			
VVT 70FG			70cm			
VVT 100FG			100cm	7cm	N/A	
VVT 70FGL*			70cm			
VVT 100FGL*			100cm			
VVT 70FGW**	7Fr	Alloy52	70cm		0.025"	
VVT 100FGW**			100cm			
VVT 4510				1cm		
VVT 4515			45cm	1.5cm		
VVT 4520				2cm		
VVT 4530				3cm		
VVT 60			60cm	7cm		
VVT 7020				2cm		
VVT 7030			70cm	3cm		
VVT 70				7cm		
VVT 90	90cm	7cm				
VVT 100	100cm	7cm				

* Led light on tip of catheter ** Embedded guidewire on tip

Compatibility of catheters with generators

		V-700	M-3004	V-1000
STS Catheter	LED	○	×	×
	General(5Fr & 7Fr) & Guide	○	×	×
Alloy Catheter	General	○	○	○